










BellMe
Saving Lives

HOME USE WIRELESS NURSE CALLING SYSTEM

031 828 2827
info@bellme.co.za
www.bellme.co.za

CONTENTS

COMPANY INFORMATION

ABOUT US.....03

HOW BELLME WORKS.....04

WHY GO WIRELESS?.....05

TESTIMONIALS..... 06

TRANSMITTER INFORMATON

B-A1-WL - SINGLE KEY CALL BUTTON.....08

RECEIVER INFORMATION

RECEIVER COMPARISON CHART.....10

B-700 LED LIGHT RECEIVER.....11

B-650 WRISTWATCH.....12

B-Q4 SIGNAL REPEATER.....13

B-99P DISPLAY PANEL.....14

CONTACT US.....15



ABOUT US

Introduction to BellMe

Founded in 2012 by Ushir Shah, BellMe came to exist in the search for a solution to simplify communication in various industries. BellMe essentially provides instant notification from cared-for to caregiver thereby assisting in reducing the reaction time in emergencies through the use of efficient and reliable Wireless Nurse Call Systems.

BellMe provides unique and revolutionary Wireless Nurse Call Systems that are fast becoming the answer to many households in South Africa.

Our Mission

To expand throughout South Africa and then Africa to become one of the largest Nurse Call System suppliers. We intend on doing this by offering our clients the best value in terms of service and product by helping you to communicate faster and more effectively, therefore reaction time in crucial situations and providing efficiency in home based care.

Our Experience

We have supplied the Akeso Group of Hospitals in South Africa, Mediclinic and Retirement Homes that span up to 600 metres. BellMe has the vision and experience to make sure we are around to take care of you and your loved ones.

OUR CREDENTIALS

ICASA APPROVED

SOLE
DISTRIBUTOR IN SOUTH AFRICA


08
AGENTS IN SOUTH AFRICA

67%
COVERAGE IN SOUTH AFRICA

12
TEAM MEMBERS

HOW BELLME WORKS

Here is a brief overview of how the system would work in a home environment:





The ***Call Button*** is mounted in high risk areas such as bedsides, showers, toilets or baths. It can also be worn on a lanyard around the neck. When pressed, a notification is wirelessly transmitted from the ***Call Button*** to:



The caregiver or a family member via a ***Display Panel*** in the lounge or any other visible area of the home. The same signal can also be sent to a ***Wristwatch*** which gives off an alert sound and a vibration in case the caregiver or family member is away from the Display Panel or during load-shedding.



The same signal may also be sent to an ***LED Light*** which can be installed in the care-giver's room or another family member's room. This light has a loud, continous alert that we recommend to be used if a care-giver or family member needs to be woken up. A ***Cancel Button*** must be present to disable the light.



Your property may need ***Signal Repeaters*** to carry the signal further or to reach areas outside the house such as a garden or garage. We will assist you in determining if you will need this product.

WHY GO WIRELESS?

Our Wireless Nurse Call Systems will significantly reduce your response time to a call for assistance in a critical situation.



Quick Installation, Easy to Use

Self Installation is a breeze. Simply mount our call buttons onto a flat surface and you are good to go.

Agents available in your area

Enquire with us for more information regarding agents in your area to conduct an on-site inspection and quotation.

No Monthly Rental

You do not need to incur extra monthly expenses when purchasing from BellMe. All our products are once-off purchases.

No Hidden Costs

Your Nurse Call System comes with the basic equipment required for it to run smoothly.

Free Support

You don't have to be slowed down by any technical aspect of our products. We will provide telephonic or remote support.

No Internet Needed

Our systems use radio frequency to wirelessly send out a signal or notification. Your internet does not affect our systems.

Solutions during load shedding

If you are constantly affected by load shedding, speak to a consultant and we will offer you a solution that works during load shedding.

The Cared-For can move freely

Wireless technology allows the cared-for to call for assistance from the bathroom, garden, lounge or any other area of the property.

Reliable Notifications

BellMe has proven to be more reliable than wired nurse call systems in the event of lightening, rodents, and electrical repairs.



TESTIMONIALS

Our experience in Nurse Call Systems speaks for itself. Have a look at what some of our clients have to say.

C

Our main reason for buying this product is because it was advertised as something that is far-reaching. We really like that the range of the product reaches far with the Signal Booster. If anyone is thinking of buying this product I suggest you buy it from Jolandie Wessels, she is very helpful and pleasant. We are very happy that the product works as advertised.

Charl Jager
MIDDELBURG, 2019

A

We are very happy with the Bellme. We got it for our boss who has a carer looking after him. He has just come home from having aggressive radium and chemo. It is working well. Thank you for your excellence service and delivery.

Ashley Hodson
PIETERMARITZBURG, 2019

R

BellMe is a wonderful product. Having parents in their late 80s or 90s is scary for any child. My father recently broke his shoulder and without the use of his right arm he found himself stuck in a couple of situations and was calling for help for over 30 minutes. Thank goodness I came across the BellMe website. I was quickly able to purchase a device for my 89 year old mother and 92 year old father. I can rest easy knowing that all they need to do is press the button on their device and their care giver is instantly notified that they need help. The good news is that BellMe is very reasonably priced though I would gladly have paid more for the peace of mind that having my parents equipped with a BellMe device has brought me.

Ralph Baard
DURBAN, 2019



TESTIMONIALS

L

My son-in-law ordered your watch/alarm for me when my wife became bed-ridden and I have found it to be most useful for keeping in touch with her when I am away from her in the house or out in the garden. It is used every day and has attracted a lot of attention from friends and potential users.

Louie van den Aardweg
DURBAN, 2019

K

Really awesome, BellMe is a really cool product, I'm actually a C6 quadriplegic and I use it to communicate with my helper, I actually have a 24-hour helper here.

Initially I use to just rely on my cell phone to contact her but you know cellphone is not a 100% solution where sometimes we have issues with the network or the phone is off and I can't always get her on time so I started using the BellMe system and I'm very very impressed with it.

Firstly it's very easy to use because I have limited finger mobility so it is easy to press the button, button is quite big. Secondly it's quite convenient, I can always have it around my neck. It is very light weighed, I don't even feel it there. I just have it around my neck the whole day and whenever I need my helper I can just call her.

The other thing is that I was concerned about was the distance because she is based quite a far distance away and it is on a level below us so I was concerned about the frequency, whether it will reach but it seems to be transmitting very well and it works 100% even though it's quite far away.

Lastly the price, I was quite impressed with the price as well, I always thought that these wireless systems were quite expensive and I found it to be quite affordable. I got 4 buttons, a calling light and a UPS. Definitely well worth it, you know it's something essential that I need in my home now especially being a quadriplegic, so it's something I can't live without.

Kreason Govender
DURBAN, 2017

B-A1-WL - SINGLE KEY CALL BUTTON



This Call Button Saves Lives.

A simple push of the button will send out alerts to multiple receivers such as:

- A Display Panel,
- A Wristwatch receiver and
- An LED Receiver.

One Button can alert a caregiver or family member through various methods.

The button is durable, waterproof and long lasting ensuring your family member has peace of mind for a lifetime.

BENEFITS



5 Year Warranty:

We have complete faith that our Buttons will continue to Save Lives.

Which is why we recently increased our Warranty on all Call Buttons from a 2-year to a 5-year warranty as of May 2018. If your Button happens to give you any issues that are covered under your warranty we will even give you a replacement unit until we resolve your technical issue.



100% Waterproof:

This call Button can be taken into high risk areas such as bathrooms, toilets and showers.

It allows the cared-for to take their trusty Call Button into the shower or out into the garden so they don't have to fear falling down and waiting for hours before someone finds them.



Dual Purpose Application:

1 - Wearable:

A lanyard can be connected to the back-mounting section and can be used around the neck by the cared-for so that they have the freedom to move around and not be confined to the area that their button is located. They can press the call button wherever they are.







A - Insert the key to release the backing.
B - The backing is mounted to a surface.
C - Call Button slides securely onto the mounted backing.

2 - Mountable

Alternatively, it can be wall mounted as it comes with a wall mounting bracket. You can drill it into the wall with the mounting bracket or use very strong 2-way tape. The Call Button can be mounted on any flat surface within the room and other high risk areas such as toilets, bathrooms and showers as they are completely waterproof.

B-A1-WL - SINGLE KEY CALL BUTTON

FEATURES

Easy to Press:

- Large surface area that can be triggered on a soft rubber with a slight press.
- The rubber surface is indented to prevent an accidental press.

Long Lasting Battery Usage:

- Requires an A23 battery.
- Battery lasts up to 1-2 years in a home-use environment.
- Supplied with a special key to remove the backing and unscrew the casing to change the battery.


Pre-programmed with a Unique Code:

- Each call button has a unique code
- The code can be saved to display what you would like.

APPLICATIONS







Showers, Bathtubs or Toilets, Bedsides:

Usually mounted in high-risk areas such as showers, bathtubs, toilets or bedsides. Can be easily reached by the cared-for and ONE simple push of the button can send an alert to a wristwatch receiver, an LED Light Receiver and the Display Panel.



Worn on a Lanyard:

If your family member also needs to move around to other areas of the house a button can also be worn on a lanyard so that they can call for assistance from any area of the house including areas such as the garden, patio or balconey or the garage. Please note, in some cases a Signal Booster may be needed however we will help you determine if your property requires a signal booster.

B-A1-WL Single Key Call Button Technical Information	Dimensions: 60x60mm	Transmission Distance: 300m in an open area
	Battery: A23	Radio Frequency: 433MHz

RECEIVER COMPARISON CHART

Please note, this is a comparison chart for quick reference to a summary of product features. It is recommended to read through the each product's benefits and features in detail to determine if it is a viable solution for your home.

Feature			
	B-99P Display Panel	**B-650 Wristwatch**	**B-700 LED Light**
Warranty	2-Years	2-Years	2-Years
Display Options	2 Numbers	1 Alphabet and 3 Numbers OR 4 Numbers	LED Light
Letters and Numbers	00-99	A001-F999 OR 0000-9999	None
Queuing Mode	No	Yes	No
Capacity	99	200	99
Continuous Ringing	No	No	Yes Must be used with Cancel Button to disable alert
Single Ring	Yes	Yes	No
Auto-Cancel	Yes[1]	No	No
Alert Tones	1	1	1
Night Use	Can be used to wake someone up at night	Not recommended	Recommended for alerts during the night
Load-Shedding	UPS needed to power on during load-shedding	Will work as long as the battery is charged	UPS needed to power on during load-shedding
Installation	Mounted on a wall OR Comes with a Display Stand	Worn on the wrist OR Around the neck with a lanyard	Mounted on a wall or above a door

[1]*You can choose between 5-99 seconds*

B-700 LED LIGHT RECEIVER

This receiver is loud enough to wake up a family member or caregiver if there is a call for assistance during the night.



BENEFITS

- Improve your response time to call for help.
- Continuous ringing until the cancel button is pressed.

FEATURES

Capacity:

This LED Light Receiver can store up to 99 call buttons & requires 1 cancel button. It can be configured by the following two methods:

- 49 Call Buttons and 49 Cancel Buttons
- 98 Call Buttons and 1 Cancel Button



Alert Tones:

Audio:

The LED has a switch inside the back cover which allows you to change between a beeping noise or silent.

Visual:

This LED light also flashes bright red upon pressing the call button.



LED Light installed above a door frame for easy identification of a Call for Assistance.

Power:

Even though it is part of our wireless range, all receivers require constant power except the wristwatch receiver. The signal transmits wirelessly from the Call Button to the LED light. The LED light is supplied with the required 12V power supply which is 1m long.

Cancel with one or multiple:

The LED light cannot work without a Cancel Button.

A cancel button is used to stop the light from flashing and making a noise. You can share one cancel button for multiple call buttons or have a cancel button for each call. The configuration would be dependent on what you like displayed on your display panel, so it is compatible with the display information required. Please speak to a consultant to find the best solution for your home.



Cancel Button must be present to disable LED Light Receiver's alert.

B-700 LED Light Receiver Technical Information	Working Voltage: 12V 500mA	Working Current: <65mA	Standby Current: <10mA	Working Environment: -40~80°C
	Sensitivity: 106dB	Volume Output: 85dB	Dimensions: 104x86x38mm	

B-650 WRISTWATCH

Do you need to receive a notifications on the move?



BENEFITS

- In a situation where you would like faster response time, not having to go to the Display Panel to view the details saves you time.
- The wristwatch gives you an immediate notification in area of the property or even if you are in the garden.
- You can cancel the request on the wristwatch, attend to the Call and then cancel at the call point.

FEATURES

Wearable on a Lanyard or Wrist:

You have two options:

Remove the watch straps and wear the wristwatch on a lanyard or wear it on your wrist.

Capacity:

The wristwatch receiver has a capacity of 200 transmitters.

The Buttons can be all Call Buttons and it has a dedicated cancel* button on the wristwatch.

**Note: Pressing the cancel button on the wristwatch will not cancel the call on the display panel. It will only clear the notification on the wristwatch display.*

Display Information:

Display either four numbers “0000-9999” or use the first digit as a letter, from A-F. Eg: “A001”-“F999”.

Rechargeable Battery:

Recharge the lithium-ion battery every 2-5 days in a nursing environment (depending on use). Battery life is indicated on the display.

Signal Strength Indication:

You can see the strength of the signal displayed on the wristwatch. If you require it to travel further than the standard 150 meters, add signal repeaters for an extended range of up to 1 km.

TECHNICAL INFORMATION



Loadshedding:

The wristwatch will work as normal during loadshedding – just ensure it is always charged. We would suggest asking your caregiver to set an alarm and stay awake during loadshedding to ensure they will be aware of any alerts.



Single Key Cancel Button:

If you use this Cancel Button (B-A1-BC) to attend to a call for assistance you will also have to cancel the call on the wristwatch using its dedicated cancel button.

B-650 Wristwatch Technical Information	Working Voltage: 3.7V (Li-ion battery)	Working Current: <900M	Frequency: 433 MHZ +/- 75K
	Working Distance: >150M	Sensitivity: 85dB	Dimensions: 65x44x19mm

B-Q4 SIGNAL REPEATER

Your property may need signal repeaters to carry the signal further.

These signal repeaters are installed where the signal falls short. It extends the range of the transmission allowing the signal to reach other areas of your property for example: the garden, a patio, a garage, etc.



BENEFITS

- Give your family member the freedom to move.
- The signal booster extends the range of the Call Button transmission.
- Your family member can have complete peace of mind knowing that their Call Button will reach you or their caregiver whether they are in the Garden, the patio or in the garage.

FEATURES

Unlimited Capacity:
The call buttons does not need to be programmed to the repeater.

Antennae:
With two Antennae to help receive and transmit the signal, the device ensures that BellMe sufficiently covers your property with ease.

Range:
Our Nurse Call Systems has a range of 100-200m without a Signal Booster. If you add Signal Boosters you can increase the range up to 1km.

Power Supply:
The Signal repeater is supplied with the required 12V adapter which is 1m long.

Power:
Even though it is part of our wireless range, all receivers require constant power except the wristwatch receiver.

Loadshedding:
In case of a power failure, the repeaters have an internal rechargeable battery pack that will power the unit for up to 1.5 – 2 hours.

TECHNICAL INFO

How do you determine if a signal repeater is needed, where to place them, and how many:
Signal repeaters are strategically placed where the signal falls short and each property is different depending on various factors. Please speak to a consultant to determine if you need one.

Extending the power supply:
Should your repeater need to be placed in a strategic position which is out of reach of any power points we can extend the Power Supply for you on request.

B-Q4 Signal Repeater Technical Information	Wireless Frequency: 433MHz	Receving & Transmitting Distance: up to 200m

B-99P DISPLAY PANEL

Ideal for Home Use.



4 easy-to-use buttons on the front.



Mounting bracket supplied to mount on a wall.



Stand supplied for table or counter-top display.

BENEFITS

- Control your Display Panel with 4 easy-to-use buttons.
- Receive alerts on a large 2-digit display which is easy to read from a distance or with a loud alert when you are not facing the display panel.

FEATURES

Display Options:
Simple 2-digit codes to locate which area the caregiver needs to go to.
The B-99P has a 2-digit display. It can display 00-99.

Stand or Mountable:
This display panel comes with a table-top stand.
You can easily place it on your counter, or if you desire it to be mounted it does also comes with a mounting bracket.

Alert Tones:
One distinct ringtone. The volume level can be adjusted to your preference.

Capacity:
Has a capacity of 99 transmitters.
This means that you can store up to 99 Call or Cancel buttons on this display.

Cancel Button Compatible:
Program cancel buttons to this display unit.
You can place cancel buttons at the call points (bedside, toilet, shower, etc) so that the caregiver can cancel the alert tone from the call point before attending to your family member.

Auto-Cancel:
Convenient and easy to cancel.
You can choose between 5-99 seconds, how long you want the calling information to stay on the screen before it cancels automatically.

Power Supply:
The B-99P is supplied with a DC12V adapter.

Loadshedding:
In the event of loadshedding, a UPS Backup Power Supply can keep the Display Panel powered on.

B-99P Display Panel Technical Information	Working Voltage: DC12V	Working Current: <450mA	Frequency: 433MHz	Working Distance: >150M
	Sensitivity: 110dB	Audio Output: 2w	Dimensions: 248x163x23mm	

CONTACT US

Please feel free to get in touch with us for any and all enquiries. We would love to hear from you.



Director
Mr. Ushir Shah
031 828 2827
ushir@bellme.co.za



Manageress
Jolandiè Wessels
031 828 2827
064 532 0828
jolandie@bellme.co.za



Admin
Terean Mooninthan
031 828 2827
064 532 0828
terean@bellme.co.za

AUTHORIZED AGENTS



Western Cape
Paul Johnson
cptsales@bellme.co.za



Kwa-Zulu Natal & Western Cape
Terence
terence@bellme.co.za



Free State, North West & Northern Cape
Martin Britz
martin@bellme.co.za



Oudtshoorn
Johan Burger
johan@bellme.co.za



Kwa-Zulu Natal
Vincent van der Walt
vincent@bellme.co.za



Eastern Cape
Ockert Davis
ockert@bellme.co.za



Mpumalanga
Barend Le Roux
barend@bellme.co.za



Gauteng
Chris Thomas
chris@bellme.co.za

Gauteng
Alan Vermaak
allan@bellme.co.za



Zimbabwe
Ngoni Wachi
ngoni@bellme.co.za

BellMe
Saving Lives

HEAD OFFICE CONTACT DETAILS

(+27) 31 828 2827
(+27) 64 532 0828
info@bellme.co.za
www.bellme.co.za
53 Herrwood Drive, Herrwood Park, Umhlanga, 4319